

Three small white hearts arranged in a slight arc above the word 'PREGNANCY'.

# PREGNANCY

RAPID TEST 



A MOMENT IN MY BELLY, A LIFETIME IN MY HEART

---

# BELLY SPOT


## PREGNANCY TEST

HCG is glycoprotein hormone secreted by the placenta shortly after implantation. In normal pregnancy, hCG can be detected in both urine and serum as early as 7 to 10 days after fertilization. The appearance of hCG and its subsequent rapid rise in urine and serum after conception during early gestational growth makes it an excellent marker.



Product name	Package	Detection limit	Specimen	Shelf life
Belly Spot HCG Test	50 Test / Box 25 Test / Box 1 Test / Box	10 mIU / ml, 25 mIU / ml	Urine / Serum	24 months

### RESULTS :



### TEST PROCEDURE :



A WOMAN WHO HAS 



A BABY FEELS A SENSE OF

**COMPLETENESS**


---



# OVULATION TEST


## LH TEST

Ovulation is the release of an egg from ovary. Luteinizing hormone (LH) which stimulates ovulation is suddenly increased (LH surge) a day before ovulation. Because the few days around ovulation are the most likely time to be pregnant, it is very important to find the time when the women ovulate. Urine LH detection is a very helpful way to conceive.




Product name	Package	Detection limit	Specimen	Shelf life
LH Card	25, 1 test / case	25 mIU / ml	Urine	24 months

### RESULTS :



### TEST PROCEDURE :





Belly  
Spot

Spot Your **BELLY**  
to Bring **HAPPINESS**

---



# MISSION

TO PROVIDE EXCELLENT  
IVD HEALTH CARE  
PRODUCT TO MEET  
HUMAN HEALTH NEEDS



RAPID TEST CATALOGUE



CONSUMABLES CATALOGUE

**SCAN FOR THE CATALOG**



INDIA | HONGKONG | TURKEY | USA |  
INDONESIA | BANGLADESH | VIETNAM



[WWW.SRBIOERA.COM](http://WWW.SRBIOERA.COM)



[INFO@SRBIOERA.COM](mailto:INFO@SRBIOERA.COM)