

SR BIOERA



— R A P I D K I T —

Rapid kit

INFECTIOUS DISEASES

PRODUCT NAME	CATALOGUE	SPECIMEN	SENSITIVITY	SPECIFICITY
Typhoid IgG/IgM Rapid Kit	SRTP01	WB/S/P	92%	99%
Dengue IgG/IgM Rapid Kit	SRDEN02	WB/S/P	97%	99.30%
Dengue NS1 Rapid Kit	SRDNS03	WB/S/P	96%	99.50%
Zika IgG/IgM Rapid Kit	SRZIK04	WB/S/P	94.10%	99.70%
CHIK IgG/IgM Rapid Kit	SRCHG05	WB/S/P	95%	99.80%
CHIK IgM Rapid Kit	SRCHM06	WB/S/P	94.30%	99.50%
Leishmania IgG/IgM Rapid Kit	SRLEG07	WB/S/P	92.40%	99.50%
Leishmania Antibody Rapid Kit	SRLEB08	WB/S/P	91.20%	99.50%
Malaria PF Antigen Rapid Kit	SRMPF09	WB/S/P	100%	100%
Malaria PF/PV Antigen Rapid Kit	SRMPV10	WB/S/P	92%	99%
Malaria Pan./PF Antigen Rapid Kit	SRMPN11	WB/S/P	93%	100%
Filaria IgG/IgM Rapid Kit	SRFGM12	WB/S/P	95.20%	100%
Filaria Antibody Rapid Kit	SRFAB13	WB/S/P	96%	100%
Leptospira IgG/IgM Rapid Kit	SRLGM14	WB/S/P	95.50%	100%
Leptospira Antibody Rapid Kit	SRLAB15	WB/S/P	96%	100%
Chagas IgG/IgM Rapid Kit	SRCGM16	WB/S/P	93%	99.60%
Chagas Antibody Rapid Kit	SRCAB17	WB/S/P	92.70%	99.20%
Malaria Pan Antigen Rapid Kit	SRMPA18	WB	-	-
Malaria Pv Antigen Rapid Kit	SRMPV19	WB	-	-
Malaria PF/PV Antibody Rapid Kit	SRMPP20	WB/S/P	-	-
Rubella IgG/IgM Rapid Kit	SRRGM21	WB/S/P	-	-
S.Typhi Antigen Rapid Kit	SRSTA22	Feces	-	-
Strep A Rapid Kit	SRSAU23	Swab	-	-
Syphilis Rapid Kit	SRSUS24	WB/S/P	-	-
Vibrio Cholerae Rapid Kit	SRVCU25	Feces	-	-

RESPIRATORY TRACT

PRODUCT NAME	CATALOGUE	SPECIMEN	SENSITIVITY	SPECIFICITY
SARS-CoV-2 IgG/IgM Rapid Kit	RSCGM01	WB/S/P	98.50%	99.90%
SARS-CoV-2 Antigen Rapid Kit	RSCB02	Dropharyngeal, Nasopharyngeal Anterior nasal (swab)	97.63%	100%
SARS-CoV-2 Antigen Rapid Kit	RSCAN03	Saliva	97.92%	98.34%
Influenza A and B Antigen Rapid Kit	RSIAB05	Dropharyngeal, Nasopharyngeal Anterior nasal (swab)	95.70%	100%
Adenovirus Antigen Rapid Kit	RSAAR06	Dropharyngeal, Nasopharyngeal Anterior nasal (swab)	94.50%	100%
Adenovirus Antigen Rapid Kit	RSAAU09	Feces	95.50%	100%
Respiratory Syncytial Virus (RSV) Antigen Rapid Kit	RSRSV07	Swab	95.50%	100%
Chlamydia pneumoniae IgM Rapid Kit	RSCP08	WB/S/P	90.50%	99%
Chlamydia pneumoniae IgG Rapid Kit	RSCP04	WB/S/P	93.20%	99.20%

SEXUALLY TRANSMITTED DISEASES

PRODUCT NAME	CATALOGUE	SPECIMEN	SENSITIVITY	SPECIFICITY
HIV				
HIV (I+ II) Antibody Test(Trilines lines) Rapid Kit	BEHIV07	WB/S/P	99.70%	99.90%
HIV (I+ II) Antibody Test(Two lines) Rapid Kit	BEHTL05	WB/S/P	99.70%	99.90%
HIV (I+II+0) Antibody Test(Trilines lines) Rapid Kit	BEHAB06	WB/S/P	99.70%	99.90%
HIV (I+II+0) Antibody Test(Two lines) Rapid Kit	BEHTA03	WB/S/P	99.70%	99.90%
HIV Ag/Ab Test Rapid Kit	BEHAG02	WB/S/P	99.70%	99.90%
HIV (I+II) Antibody Test(Two lines) Rapid Kit	BEHTD4	SALIVA		
HIV (I+ II) Antibody Test(Two lines) Rapid Kit	BEATB01	URINE		
TP				
TP Antibody Rapid Kit	BETPA10	WB/S/P	100%	98.50%
HERPES SIMPLEX VIRUS-I/II (HSV-I/IL)				
HSV-I IgM Rapid Kit	BEHSI16	WB/S/P	91.20%	99%
HSV-I IgG Rapid Kit	BEHSL13	WB/S/P	94.20%	99.50%
HSV-I IgG /IgM Rapid Kit	BEHGM12	WB/S/P	93.60%	99%
HSV-II IgM Rapid Kit	BEHIM15	WB/S/P	90.20%	99.10%
HSV-II IgG Rapid Kit	BEIG11	WB/S/P	91.20%	99%
HSV-II IgG/IgM Rapid Kit	BEHGG14	WB/S/P	93.60%	99%
COMBO TEST				
HIV/HCV (trilines) Antibody Rapid Kit(Trilines lines)	BEHVA18	WB/S/P	99.70%	99.80%
HIV/TP(four-lines) Antibody Rapid Kit(Trilines lines)	BEHTP17	WB/S/P	99.70%	99.50%

CARDIAC MARKERS

PRODUCT NAME	CATALOGUE	SPECIMEN	SENSITIVITY
Cardiac Troponin I Rapid Kit	SBCDM02	WB/S/p	0.5 ng/ml
CK-MB Rapid Kit	SBCDM01	S/p	5 ng/ml

VETERINARY DIAGNOSTICS


PRODUCT NAME	CATALOGUE	SPECIMEN
CANINE DISEASE RAPID KIT		
CPV Antigen Rapid Kit	ARCPV14	Body Secretion
Canine CRP Antigen Rapid Kit	ARCR10	Feces
CDV Antigen Rapid Kit	ARCDV13	Feces
CCV Antigen Rapid Kit	ARCCV12	Body Secretion
Canine Influa Antigen Rapid Kit	ARCIA09	Feces
CCV/CPV Antigen Rapid Kit	ARCPV11	Body Secretion
Rota virus Antigen Rapid Kit	ARRAV08	
CHW Antigen Rapid Kit	ARCHW07	WB/S/P
CDV Antibody Rapid Kit	ARCDV05	WB/S/P
CPV Antibody Rapid Kit	ARCPV06	WB/S/P
CPV/CDV Antibody Rapid Kit	ARCAB01	WB/S/P
FELINE DISEASE RAPID KIT		
FPV Antigen Rapid Kit	ARFPV03	Feces
FIV Antigen Rapid Kit	ARFIV02	WB/S/P
FelV Antigen Rapid Kit	ARFLV04	WB/S/P
PARASITIC DISEASE AND OTHERS		
TOXD Antibody Rapid Kit	ARTAB15	WB/S/P
Brucella Antibody Rapid Kit	AR	WB/S/P
Rota virus/Adenovirus Antigen U.S.	AR	Feces

DRUGS OF ABUSE

PRODUCT NAME	CATALOGUE	SPECIMEN	SENSITIVITY
ACE Acetaminophen Rapid Kit	SS00A01	URINE/SALIVA	5000 ng/ml
AMP Amphetamine Rapid Kit	SS00A02	URINE/SALIVA	500/1000 ng/ml
BAR Barbiturates Rapid Kit	SS00A03	URINE/SALIVA	300 ng/ml
BUP Buprenorphine Rapid Kit	SS00A04	URINE/SALIVA	10 ng/ml
BZO Benzodiazepines Rapid Kit	SS00A05	URINE/SALIVA	200/300 ng/ml
CAF Caffeine Rapid Kit	SS00A06	URINE/SALIVA	10000 ng/ml
COC Cocaine Rapid Kit	SS00A07	URINE/SALIVA	100/200/300 ng/ml
COD Codeine Rapid Kit	SS00A08	URINE/SALIVA	300 ng/ml
COT Cotinine Rapid Kit	SS00A09	URINE/SALIVA	100/200/500 ng/ml
EDDP Methadone Metabolite Rapid Kit	SS00A10	URINE/SALIVA	100 ng/ml
ETG Ethyl Glucuronide Rapid Kit	SS00A11	URINE/SALIVA	1000 ng/ml
FYL Fentanyl Rapid Kit	SS00A12	URINE/SALIVA	200 ng/ml
KET Ketamine Rapid Kit	SS00A13	URINE/SALIVA	500/1000 ng/ml
Synthetic Marijuana (K2) Rapid Kit	SS00A14	URINE/SALIVA	50/200 ng/ml
MCA Methcathinone Rapid Kit	SS00A15	URINE/SALIVA	1000 ng/ml
MDMA (Ecstasy) Rapid Kit	SS00A16	URINE/SALIVA	500 ng/ml
MDPV Rapid Kit	SS00A17	URINE/SALIVA	1000 ng/ml
MET Methamphetamine Rapid Kit	SS00A18	URINE/SALIVA	300/500/1000 ng/ml
MOP Morphine Rapid Kit	SS00A19	URINE/SALIVA	100/200/300 ng/ml
MLQ Methaqualone Rapid Kit	SS00A20	URINE/SALIVA	300 ng/ml
MTD Methadone Rapid Kit	SS00A21	URINE/SALIVA	300 ng/ml
OPI Opiate Rapid Kit	SS00A22	URINE/SALIVA	300/2000 ng/ml
OXY Oxycodone Rapid Kit	SS00A23	URINE/SALIVA	100 ng/ml
PCP Phencyclidine Rapid Kit	SS00A24	URINE/SALIVA	25 ng/ml
PPX Propoxyphene Rapid Kit	SS00A25	URINE/SALIVA	300 ng/ml
Tricyclic Antidepressants (TCA) Rapid Kit	SS00A26	URINE/SALIVA	1000 ng/ml
Marijuana (THC) Rapid Kit	SS00A27	URINE/SALIVA	25/50 ng/ml
TRA Tramadol Rapid Kit	SS00A28	URINE/SALIVA	200/300 ng/ml
Alcohol Rapid Kit	SS00A29	URINE/SALIVA	See Insert

GASTROINTESTINAL PATH PGENS

PRODUCT NAME	CATALOGUE	SPECIMEN	SENSITIVITY	SPECIFICITY
H.Pylori Antigen Rapid Kit	MDHPA03	WB/S/P	94.40%	100%
H.Pylori Antibody Rapid Kit	MDHAB01	WB/S/P	90%	95%
F08 Antigen Rapid Kit	MDFOB02	WB/S/P	50 ng hHb/mL	-





TORCH & CHILDHOOD

PRODUCT NAME	CATALOGUE	SPECIMEN	SENSITIVITY	SPECIFICITY
TOXD IgM Rapid Kit	RRTOX05	WB/S/P	91.60%	99%
TOXD IgG Rapid Kit	RRIG607	WB/S/P	93%	99.20%
TOXD IgG/IgM Rapid Kit	RRIGM04	WB/S/P	91.80%	99%
CMV IgM Rapid Kit	RRCGM01	WB/S/P	92.70%	99.10%
CMV IgG Rapid Kit	RRCGG02	WB/S/P	93.60%	99%
CMV IgG/IgM Rapid Kit	RRCGM06	WB/S/P	93%	99.20%
Rota virus Antigen Rapid Kit	RRVRA03	Swab,Feces	99.70%	99.90%
HAV IgG/IgM Rapid Kit	RRHAV08	WB/S/P	-	-
HBcAb Rapid Kit	RRHBC09	WB/S/P	-	-
HBeAg Rapid Kit	RRHBE10	WB/S/P	-	-
HCV/HBsAg/HIV Combo Rapid Kit	RRHHH12	WB/S/P	-	-
HBV Rapid Kit(HBsAg/HBsAb/HBeAg/HBeAb/HBcAb)	RRHBV13	WB/S/P	-	-
HEV IgG/IgM Rapid Kit	RRHEV14	WB/S/P	-	-

TUMOR MARKERS

PRODUCT NAME	CATALOGUE	SPECIMEN	SENSITIVITY
Fecal Occult Blood Rapid Kit	SCTOM01	Feces	10/40/50/100/200 ng/ml

HEPATITIS


PRODUCT NAME	CATALOGUE	SPECIMEN	SENSITIVITY	SPECIFICITY
HCV Antibody Rapid Kit	SDHCV01	WB/S/P	100%	99.20%
HBs Antigen Rapid Kit	SDHBS02	WB/S/P	96%	100%
HAV IgM Rapid Kit	SDHAV03	WB/S/P	96%	99.20%
HEV IgM Rapid Kit	SDHEV04	WB/S/P	99.70%	99.90%
HBsAg/HCV (trilines) Rapid Kit	SOSA605	WB/S/P	96%	99.20%

FERTILITY

PRODUCT NAME	CATALOGUE	SPECIMEN	SENSITIVITY
HCG Pregnancy Rapid Kit	SZFER05	URINE	25 mIU/ml
HCG Pregnancy Rapid Kit	SZFER03	URINE/SERUM	10/25 mIU/ml
LH Ovulation Rapid Kit	SZFER01	URINE	25 mIU/ml
FSH Menopause Rapid Kit	SZFER02	URINE	25 mIU/ml
Sperm SP10 Protein Rapid Kit	SZFER04	SEMEN	15x10 ⁷ ml

A COMPANY OF
QUALITY

A PLACE
WITH A MOTIVE



DELIVERING

AFFORDABLE

INNOVATION

MISSION



TO PROVIDE EXCELLENT
IVD HEALTH CARE
PRODUCT TO MEET
HUMAN HEALTH NEEDS



RAPID TEST CATALOGUE



CONSUMABLES CATALOGUE

SCAN FOR THE CATALOG



PL E-30, HINDWA DREAMS,
DHORAN-PARDI, KAMREJ-
394155 SURAT, GUJARAT, INDIA.



WWW.SRBIOERA.COM



INFO@SRBIOERA.COM